


Outline

- LEAP results on Nitronex devices
 - Characterization of gate/epilayer interfacial layer
- Challenges with performing LEAP
 - LEAP results on AFRL devices
 - Sample fracture during LEAP
 - Interfacial layer thickness variations
- Addressing LEAP challenges
 - Adjustment of laser wavelength
- Update on collaboration
 - Dr. Hono's group (University of Tsukuba, Japan)
- Conclusions and future work

TEM of Gate/Epi-layer Interface (Nitronex)


HAADF-STEM Images




- Interfacial layer ~ 0.5 nm-thick evident between Ni and AlGaN; can perform LEAP to determine composition

LEAP of Gate Stack and Epilayers (Nitronex)

HAADF-STEM




LEAP




- Gate metal stack, gate/epilayer interface, AlGaIn/GaN interface contained in LEAP reconstruction

LEAP Reconstruction Movie (Nitronex)




Analysis Cross Section (Nitronex)

LEAP reconstruction showing all atoms



- Position a “data pipe” orthogonal to the layers to produce 1D concentration profiles for each species


1D Interface Concentration Profiles (Nitronex)



- Perform closer examination of the Ni/AlGaN interface to characterize interfacial layer

M.R. Holzworth, *et al.*, Appl. Phys. Lett. **98**, 122103 (2011).


A Closer Inspection of the Interface (Nitronex)



- Interfacial layer composed of distinct NiO_x and AlO_x layers, but what about GaO_x ?

M.R. Holzworth, *et al.*, Appl. Phys. Lett. **98**, 122103 (2011).


Is GaO_x Present in the Interfacial Layer? (Nitronex)



- Expect peaks at 85.0 and 87.0 for GaO , but none observed; should also look at 42.5 and 43.5 for doubly-ionized species

M.R. Holzworth, *et al.*, Appl. Phys. Lett. **98**, 122103 (2011).

Is GaO_x Present in the Interfacial Layer? (Nitronex)



- $^{69}\text{GaN}^{++}/^{71}\text{GaN}^{++} = 1.52 \pm 0.05$ (close to $^{69}\text{Ga}/^{71}\text{Ga} = 1.51$) and lack of peak at 43.5 indicate no GaO_x present



M.R. Holzworth, *et al.*, Appl. Phys. Lett. **98**, 122103 (2011).

Outline

- LEAP results on Nitronex devices
 - Characterization of gate/epilayer interfacial layer
- Challenges with performing LEAP
 - LEAP results on AFRL devices
 - Sample fracture during LEAP
 - Interfacial layer thickness variations
- Addressing LEAP challenges
 - Adjustment of laser wavelength
- Update on collaboration
 - Dr. Hono's group (University of Tsukuba, Japan)
- Conclusions and future work

TEM of Gate/Epi-layer Interface (AFRL)


HAADF-STEM Images




- Interfacial layer ~1.5 nm-thick present in AFRL devices; attempted LEAP analysis to determine composition

LEAP of AFRL Devices

HAADF-STEM




LEAP



- Ni and AlGaN layers thinner in LEAP than observed in TEM; indicates fracture at Ni/AlGaN interface during LEAP

Mass/Charge History Comparison




- Nitronex: M/C history broad and graded; indicates continuous evaporation without fracture
- AFRL: M/C history shows sharp discontinuity at interface; indicates tip fracture at gate/epilayer interface

Limitations of LEAP for Analyzing Devices

- Fracture during LEAP analysis of ARFL and TriQuint devices occurs near interfaces
- Possible explanations for fracture at interfaces
 - Weak interfacial adhesion
 - Interfacial layer thickness
 - Stresses generated during laser pulsing
 - Gate morphology (cracks)
- LEAP was successful on Nitronex, but not AFRL or TriQuint devices, so what is inherently different?


Differences in Interfacial Layers Between Devices

Nitronex




~0.5 nm

AFRL



~1.5 nm

TriQuint





~1.8 nm

- Nitronex HEMTs have smallest interfacial layer; these were the only devices that didn't fracture during LEAP analysis

Gate Morphology Effects: Cracks (AFRL)

HAADF-STEM Images



- Current atom probe samples employ top-down serial sampling method: Au \rightarrow Ni \rightarrow AlGaIn \rightarrow GaN
- Cracks in the gate contacts could cause fracture during LEAP analysis (just a speculation)

Outline

- LEAP results on Nitronex devices
 - Characterization of gate/epilayer interfacial layer
- Challenges with performing LEAP
 - LEAP results on AFRL devices
 - Sample fracture during LEAP
 - Interfacial layer thickness variations
- Addressing LEAP challenges
 - Adjustment of laser wavelength
- Update on collaboration
 - Dr. Hono's group (University of Tsukuba, Japan)
- Conclusions and future work

Addressing LEAP Challenges


- Use APT systems with shorter wavelength lasers
 - LEAP limited to laser wavelength of 532 nm
 - Newer systems have wavelengths < 360 nm
 - Small wavelength = higher photon energy
- Small wavelength improves evaporation during APT
 - Direct excitation (bandgap $<$ photon energy)
 - Thermal excitation (bandgap $>$ photon energy)
 - Both result in enhanced tip conductivity, improved evaporation
- Recently established collaboration with short wavelength laser APT capabilities
 - Dr. Hono (University of Tsukuba, NIMS at MMU)

Atom Probe Collaboration


- Dr. Kazuhiro Hono
 - University of Tsukuba (Tsukuba, Japan)
 - NIMS & MMU (Magnetic Materials Unit)
- Custom built atom probe system
 - Four laser wavelengths: 1030, 515, 343, 258 nm
 - Wavelength for present work limited to 532 nm
 - Demonstrated APT results on insulating/ wide band gap materials
 - ZnO, MgO, LiCoNiMgAlO₂ (Li-ion battery)
 - InGaN/GaN Quantum Wells
- Ray currently in Japan working w/ Hono's group

Some Results from Hono's Group

HAADF-STEM



APT Reconstruction



- Hono's APT system capable of analyzing InGaN/GaN multi-quantum well structures (similar to HEMT materials)

S. Tomiya *et al.*, Appl. Phys. Lett. **98**, 181904 (2011).

Conclusions

- 0.5 nm-thick interfacial layer in Nitronex devices analyzed using LEAP
 - Identified as distinct NiO_x and AlO_x layers
 - GaO_x not present
- Difficulties analyzing AFRL and TriQuint devices
 - Thicker interfacial layers compared to Nitronex devices
 - Cracks in gates of AFRL devices
- Started collaboration with Hono's group
 - System may be effective with thicker interfacial layers
 - Samples must be made on-site for best results

Future Work

- Atom probe collaboration with Hono's group
 - Currently in progress
 - Analyze structure/composition of defects in AFRL devices
- Structure-composition relationships of defects
 - Provide insights into formation mechanisms
- Correlation of electrical degradation with defect structure, composition, etc.